

pRP^{iTAG}

For inducible expression of N/C-terminal human cMYC- (EQKLISEEDL) or N/C-terminal enhanced GFP-tagged proteins in *T. brucei* from a tetracycline-responsive RRNA promoter.

- High fidelity polymerase recommended.
- Should integrate into any *T. brucei* genome (*RRNA*) following digestion with *NotI*.

Cloning

GFP/MYC_X

To ensure that your gene is in frame with the tag, place the second codon downstream of the *XbaI* site, i.e. TCTAGA:[codon 2]:[codon]_n:[stop codon]:GGATCC.

X^{GFP}

To ensure that your gene is in frame with the tag, place the last but one codon upstream of the *XbaI* site, i.e. AAGCTT:[start]:[codon]_{n-1}:TCTAGA.

X^{6MYC}_X

N-terminal tagging: clone your gene without a start codon via *AvrII/BamHI* digestion (or without start/stop codons via *AvrII*)
i.e. CCTAGG[codon 2]:[codon]_n:[stop]:GGATCC.

C-terminal tagging: clone your gene without a stop codon via *HindIII*(or *PacI*)/*XbaI* digestion (without start/stop codons via *XbaI*)
i.e. AAGCTT:[start]:[codon]_n:TCTAGA.

Use *HindIII* / *BamHI* if you don't want the tag.

There are alternatives if the gene contains *XbaI*, *AvrII* and/or *BamHI*:

Plasmid	<i>XbaI</i> , <i>AvrII</i> <i>BamHI</i>	Insert	<i>AvrII</i> , <i>NheI</i> , <i>SpeI</i> , <i>XbaI</i> <i>BglII</i>
----------------	--------------------------------------------	---------------	------------------------------------------------------------------------

Key features

- Complete sequences available.
- Hygromycin (GFP_X, MYC_X, X^{6MYC}_X) or Blasticidin (MYC_X, X^{6MYC}_X, X^{GFP}) for stable selection.
- All vectors allow inducible expression using tetracycline (or analogues).
- Inducible cassette is independent of selectable marker.
- Modular nature allows tag or other components to be exchanged.
- Compatible with wild type cells (for constitutive expression) and any *T. brucei* cell line expressing TetR.

Upon integration into *T. brucei*, the construct lies between *RRNA* spacer sequences. The operator binds Tet-repressor in the absence of tetracycline so the inducible RRNA promoter is activated and tagged protein is expressed in the presence of tetracycline (1 µg ml⁻¹).

New technology development

- These vectors transform bloodstream-form cells at a low efficiency (~10⁻⁷) and can integrate in the genome at any one of several ribosomal spacer loci. Subsequent position effects can generate variable results. Integrating at a tagged locus can alleviate these effects and improve transformation efficiency (see pRPa and Alsford et al, 2005).

Detection:

cMYC	Mouse anti-cMYC, 9E-10 (Source Biosciences; IFA / western blotting) Mouse anti-cMYC, 4A6 (Upstate Biotech; WB only; in IFA binds <i>Tb</i> spindle)
eGFP	Rabbit anti-GFP, (Molecular Probes; IFA and western blotting)

Other questions/comments, contact Sam Alsford (sam.alsford@lshtm.ac.uk).

pRP^{HYG-iSL}

For inducible expression of stem loop RNA in *T. brucei* from a tetracycline-responsive RRNA promoter.

- Primers: We use a software tool (RNAit – Redmond et al, 2003) for the selection of RNAi targets that provides primer information and minimises off-target effects (see Durand-Dubief et al, 2003).
- Integrates at a random RRNA spacer in *T. brucei* after *NotI* digestion.


Primer Design & Cloning

A single primer pair is required to generate the two inserts. Each primer contains two restriction sites; the internal sites allow cloning in a sense orientation into MCS1 or 3, while the external sites allow cloning in the antisense orientation into MCS2 or 4.

For example (cloning into pRP^{HYG-iSL} MCS^{1/2}):

Primer A: 5'-GATC GGGCCC GGTACC -- target specific 5' sequence (20 bases) --
Bsp120I *Acc65I*

Primer B: 5'-GATC TCTAGA GGATCC -- target specific 3' sequence (20 bases) --
XbaI *BamHI*



To confirm the organisation of the stem loop cassette, use the sense fragment cloning restriction enzymes. In the above example, *Acc65I* will excise both fragments and *lacZ*⁴⁶⁸, while *BamHI* will only release *lacZ*⁴⁶⁸.

We regularly use two sequencing primers to confirm correct insertion:

<i>lacZ</i> ⁴⁶⁸ pos 50 towards sense RNAi	Seq1 5'-AATAGTGGACTCTTGTTC
<i>lacZ</i> ⁴⁶⁸ pos 420 towards antisense RNAi	Seq2 5'-AAAGGGGATGTGCTGCAAG

Key features

- Complete sequences available.
- Hygromycin for stable selection.
- All vectors allow inducible expression using tetracycline (or analogues).
- Inducible cassette is independent of selectable marker.
- Modular nature allows components to be exchanged.
- Compatible with any *T. brucei* cell line expressing TetR.
- Upon integration into *T. brucei*, the construct lies between RRNA spacer sequences. The operator binds Tet-repressor in the absence of tetracycline so the inducible RRNA promoter is activated and tagged protein is expressed in the presence of tetracycline (1 µg ml⁻¹).

New technology development

- These vectors transform bloodstream-form cells at a low efficiency (~10⁻⁷) and can integrate in the genome at any one of several ribosomal spacer loci. Subsequent position effects can generate variable results. Integrating at a tagged locus can alleviate these effects and improve transformation efficiency (see pRPa and Alsford et al, 2005).

Other questions/comments, contact Sam Alsford (sam.alsford@lshtm.ac.uk).