





Laboratory Working Practice Document 3: Sample Processing and Storage

Ambition Trial Coordinating Centre
Private Bag 320, Princess Marina Hospital
Gaborone, Botswana

Laboratory Work Practice Document: 3 (Cape Town) Sample Processing and Storage

| | | | |
|----------------------------|--|---|------------|
| Title of study | High Dose AMBISOME® on a Fluconazole Backbone for Cryptococcal Meningitis Induction Therapy in sub-Saharan Africa: A Phase III Randomized Controlled Non-inferiority Trial | | |
| Acronym | Ambition-cm – AMBIsome Therapy Induction Optimization | | |
| ISRCTN No.: | ISRCTN72509687 | | |
| WPD Current version | Version 1.2 13/09/2019 | | |
| Author(s) | Kwana Lechiile Lab Scientist |  | 13/09/2019 |
| | Timothée Boyer Chammard Clinical Advisor |  | 13/09/2019 |
| Reviewer(s) | David Lawrence Lead Clinician |  | 13/09/2019 |
| Approved by | Joseph Jarvis CI |  | 13/09/2019 |

| Revision History: | | |
|-------------------|----------------|--|
| Version Number | Effective Date | Reason for Change |
| 1.0 | 14/06/2018 | First version |
| 1.1 | 05/12/2018 | Mbarara added to the sites, Harare removed |
| 1.2 | 13/09/2019 | <ul style="list-style-type: none"> - Lilongwe, Kampala and Mbarara removed - Addition of the qPCR sub-study in Cape Town |
| | | |
| | | |

Laboratory Working Practice Document 3: Sample Processing and Storage

Purpose

This document describes the day to day laboratory duties in Cape Town

References

Ambition Phase III Trial Protocol

Materials

Lab WPD 1 – Laboratory equipment
Lab WPD 2 – Quantitative Cryptococcal cultures
Lab WPD 7 – Quantitative Cryptococcal PCR sub-study
Clinical WPD 15- Timing of evaluation and tests
Clinical WPD 20 – Quantitative PCR Sub-study
Cape Town qPCR sub-study protocol

Appendices

Lab Processing Chart: Cape Town

Section A – Lab maintenance, quality control

Monitoring of freezer and fridge temperatures, reagents/agar expiration dates, working order of autoclaves, incubators and other laboratory equipment, should all be performed in line with on-site laboratory quality control and health and safety procedures

Weekly stock check of lab consumables for trial. Inform trial team on-site and trial manager when replacements required. See Laboratory WPD 1 – Laboratory equipment

Laboratory Working Practice Document 3: Sample Processing and Storage

Section B – Samples required

All results should be recorded with patient study number, date, and study day.

All results should be communicated to the study team.

– Study bloods

The following must be performed according to the trial protocol

- CD4 at baseline if not done within 3 months
- FBC and ALT on days 1, 7, 14, 28
- Creatinine, electrolytes on days 1, 3, 5, 7, 10, 12, 14, 28
- qPCR samples (4ml of blood – EDTA tubes): D1 (pre-dose), D3, D7 and D14,
(please see Appendix, Lab WPD7 & Clinical WPD 20 - Quantitative PCR Sub-study)

– CSF

- India Ink examination and Cryptococcal antigen (IMMY CrAg) on D1
- CSF cell count and differential, protein, glucose and routine culture on D1 for all patients and additional LPs, as indicated
- Quantitative fungal culture will be conducted on D1, 7, and 14
- CSF for qPCR (1ml): D1 (pre-dose), D7 and D14
(please see Appendix, Lab WPD7 & Clinical WPD 20 - Quantitative PCR Sub-study)

– Sputum

- Sputum for qPCR at D1 or D2

– Urine

- Urine for qPCR (20 -50 ml): D1, D3 and D7

Laboratory Working Practice Document 3: Sample Processing and Storage

Table 1: Timings of study bloods and CSF analyses

| Event Schedule | Screening | Week 1 | | | | | | | Week 2 | | | | | | Wk 4 | |
|-------------------------------------|-----------|--------|----|----|----|----|----|----|--------|----|-----|-----|-----|-----|------|---|
| Study Day | ≤D0 | D1 | D2 | D3 | D4 | D5 | D6 | D7 | D8 | D9 | D10 | D11 | D12 | D13 | D14 | |
| Clinical labs | | | | | | | | | | | | | | | | |
| HIV testing (If status unknown)* | X | | | | | | | | | | | | | | | |
| Pregnancy Test (Urine/Serum) (1) | X | | | | | | | | | | | | | | | |
| Full Blood Count | X | | | | | | | X | | | | | | | X | X |
| CD4 count & Viral load (If needed)* | X | | | | | | | | | | | | | | | |
| ALT | X | | | | | | | X | | | | | | | X | X |
| Urea, creatinine and electrolytes | X | | X | | | X | | X | | X | | X | | X | X | X |
| Blood for drug levels (2) | X | | | | | | | | X | | | | | | | |
| CSF | | | | | | | | | | | | | | | | |
| Opening pressure | X | | | | | | | X | | | | | | | X | |
| Cell count and differential* | X | | | | | | | | | | | | | | | |
| Protein, glucose* | X | | | | | | | | | | | | | | | |
| Routine culture * | X | | | | | | | | | | | | | | | |
| India Ink examination* (3) | X | | | | | | | | | | | | | | | |
| Cryptococcal antigen* (3) | X | | | | | | | | | | | | | | | |
| Quantitative fungal culture | X | | | | | | | X | | | | | | | X | |
| CSF Drug levels (2) | X | | | | | | | X | | | | | | | X | |
| Immune parameters (2) | X | | | | | | | X | | | | | | | X | |

*Part of routine care. 1. For women of childbearing age. 2. Sub-studies at limited sites 3. India ink or CrAg required for inclusion.

Laboratory Working Practice Document 3: Sample Processing and Storage

Section C - Procedures for sample processing and saving - Please see Appendix: Lab Processing Chart

Day 1 (or 0)

Required samples:

- 2 yellow (or red) top (1 for biochemistry, 1 for serum storage)
- ≥3 purple top (1 for FBC, 1 for whole blood storage, 1 for plasma and buffy coat storage and other samples may be for CD4, viral load and resistance testing)
- CSF sample in 3 x white top tube (5 ml each)
- qPCR samples

(Please see WPD 20 – Quantitative PCR sub-study)

DAY 1 CSF SAMPLES

NB CSF must be processed within **4 HOURS** of receipt

1. Perform protein and glucose on CSF
2. Perform India Ink, CrAg (if required), cell count (on day 0/1 sample)
3. Plate out day 1 unspun CSF dilutions 0-4 (neat, one in 10, 100, 1,000 and 10,000) for quantitative culture (see WPD 2)

Store isolate once grown in 2 Microbank vials.
4. Vortex CSF and save: divide into 3 x 2ml Cryovials - min 200 µL each. Label and place into freezer box.
5. qPCR SAMPLE: Please see Lab WPD7 & Clinical WPD 20 – Quantitative PCR Sub-study

TB culture only if clinically indicated (will be performed in TB lab, and is not mandated by trial protocol)

DAY 1 BLOOD SAMPLES

1. Take blood for CD4, viral load and/or resistance testing (if indicated), whole blood, plasma and buffy coat (all purple tops) and serum (yellow or red top). Total required = 2 yellow or red top (4 ml) and ≥3 purple tops (4 ml each)
2. Perform FBC, ALT, Creatinine and electrolytes on Day 1
3. **WHOLE BLOOD:** take 2x cryovials and pipette 900 µL of WHOLE BLOOD (unspun sample) into each (should be mixed well in purple top tubes prior to transfer). Label and place into whole blood storage box.
4. **PLASMA AND BUFFY COAT:** Spin blood from purple top tube at 800xg for 15 minutes. Avoiding the buffy coat and red cells at the bottom, remove as much of the plasma as possible using a sterile Pasteur pipette and divide into 2 cryovials.

Buffy coat is the small region between the plasma and red blood cells. After removing the plasma, use a sterile Pasteur pipette to carefully aspirate the buffy coat in a circular movement trying not to disturb the red blood cell layer then store in a single cryovial

Laboratory Working Practice Document 3: Sample Processing and Storage

5. **SERUM:** Spin yellow (or red) top sample at 800xg for 15 minutes. In the yellow top, gel now separates cells from serum. . Use a sterile Pasteur pipette to remove as much of the serum possible and divide into 2 cryovials.
6. **qPCR SAMPLE:** Please see Lab WPD7 & Clinical WPD 20 – Quantitative PCR Sub-study

DAY 1 SPUTUM SAMPLE

1. Process with dithiothreitol according to manufacturer's instructions (e.g. volume/volume and incubation 10min at 37°C).
2. Centrifugation at 10,000 g for 3-5 min
3. Remove the supernatant and discard
4. Transfer the pellet into an 1-2ml Eppendorf tube
5. Store at -80°C

DAY 1 URINE SAMPLE

1. Centrifuge urine sample at 4000 rpm
2. Remove the supernatant and discard
3. Transfer the pellet into an 1-2ml Eppendorf tube
4. Store at -80°C

Day 3

Required samples:

- 1 purple top (for plasma and buffy coat)
- 1 yellow/red top tube (4ml)
- qPCR samples

(Please see Lab WPD7 & Clinical WPD 20 – Quantitative PCR sub-study)

DAY 3 BLOOD SAMPLES

1. Perform urea, creatinine and electrolytes on Day 3
2. PLASMA AND BUFFY COAT
3. **qPCR SAMPLE:** Please see Lab WPD7 & Clinical WPD 20 – Quantitative PCR Sub-study

Laboratory Working Practice Document 3: Sample Processing and Storage

DAY 3 URINE SAMPLE

1. Centrifuge urine sample at 4000 rpm
2. Remove the supernatant and discard
3. Transfer the pellet into an 1-2ml Eppendorf tube
4. Store at -80°C

Day 5

Required samples:

- 1 yellow/red top tube (4ml)

DAY 5 BLOOD SAMPLES

1. Perform urea, creatinine and electrolytes on Day 5

Day 7

Required samples:

- 2 purple top (1 for FBC, 1 for plasma and buffy coat storage)
- 1 yellow/red top tube (4ml)
- CSF sample in 2 x white top tubes (5 ml each)
- qPCR samples (Please see Lab WPD7 & Clinical WPD 20 – Quantitative PCR sub-study)
-

DAY 7 CSF SAMPLES

1. Perform quantitative cultures on Day 7 LP.
2. Save remaining CSF, divide into 3 x 2ml Cryovials.
3. Save Cryptococcus isolate in 2 Microbank vials if growth from Day 7 LP
4. qPCR SAMPLE: Please see Lab WPD7 & Clinical WPD 20 – Quantitative PCR Sub-study

DAY 7 BLOOD SAMPLES

1. Perform FBC, ALT, urea, creatinine and electrolytes on Day 7
2. PLASMA AND BUFFY COAT
3. qPCR SAMPLE: Please see Lab WPD7 & Clinical WPD 20 – Quantitative PCR Sub-study

Laboratory Working Practice Document 3: Sample Processing and Storage

DAY 7 URINE SAMPLE

1. Centrifuge urine sample at 4000 rpm
2. Remove the supernatant and discard
3. Transfer the pellet into an 1-2ml Eppendorf tube
4. Store at -80°C

Day 10

Required samples:

- 1 yellow/red top tube (4ml)

DAY 10 BLOOD SAMPLES

1. Perform urea, creatinine and electrolytes on Day 10

Day 12

Required samples:

- 1 yellow/red top tube (4ml)

DAY 12 BLOOD SAMPLES

1. Perform urea, creatinine and electrolytes on Day 12

Day 14

Required samples:

- 2 purple top (1 for FBC and 1 for plasma and buffy coat)
 - 1 yellow/red top tube (4ml)
 - CSF sample in 2 x white top tube (5 mls each)
 - qPCR samples (including India ink storage)
- (Please see Lab WPD7 & Clinical WPD 20 – Quantitative PCR sub-study)

DAY 14 CSF SAMPLES

1. Perform quantitative cultures on Day 14 LP. Save remaining CSF, divide into 3 x 2ml Cryovials.
2. Save Cryptococcus isolate in 2 Microbank vial if growth from Day 14 LP
3. qPCR SAMPLE: Please see WPD 20 – Quantitative PCR Sub-study

Laboratory Working Practice Document 3: Sample Processing and Storage

DAY 14 BLOOD SAMPLES

1. Perform FBC, ALT, urea, creatinine and electrolytes on Day 14
2. PLASMA AND BUFFY COAT
3. qPCR SAMPLE: Please see WPD 20 – Quantitative PCR Sub-study

Day 28

Required samples:

- 1 purple top (for FBC)
- 1 yellow/red top tube (4ml)

DAY 28 BLOOD SAMPLES

1. Perform FBC, ALT, urea, creatinine and electrolytes on Day 28
- 2.

PROCESSING AND LABELLING

SPINNING: For all samples use 800xg for 15 mins, unless stated otherwise on LPC.

LABELLING: Each sample must be labelled with study number (e.g. AMB 1-1-001), nature of sample (e.g. CSF/plasma etc), study day (e.g. d3), date and time (time specimen taken, not time specimen received in lab).

The study day is derived from **the first day the patient receives study drugs**. That is day 1. That is not necessarily the day of hospital admission.

STORAGE: Cryovials containing samples should be stored in appropriately labelled cryoboxes. Cryoboxes should be stored in -80°C freezer

CRYPTOCOCCAL ISOLATE STORAGE: Cryptococcal isolates are stored on Microbank beads in an appropriately labelled cryobox and stored in -80°C freezer.

PROCEDURE FOR STORING ON MICROBANK BEADS:

1. Using a permanent marker, label with the study number, date and day of study (e.g. AMB 1-1-001, 12.10.15, d7)
2. Under aseptic conditions open the screw cap cryovial.
3. Inoculate the cryopreservative fluid with young colonial growth (18-24 hours) picked from a pure culture to approximately a 3-4 McFarland standard.
4. Close vial tightly and invert 4-5 times to emulsify organism. **DO NOT VORTEX!**
5. At this point the microorganisms will be bound to the porous beads. The excess cryopreservative fluid should be well aspirated leaving the inoculated beads as free of liquid as possible. Close the vial finger tight.
6. Record on the storage log.
7. Store the inoculated cryovial in -80°C freezer.

Laboratory Working Practice Document 3: Sample Processing and Storage

Section D - Summary of samples to be saved

- Day 1:**
- Serum (stored in 2 Cryovials)
 - Plasma (divided into 2 Cryovials)
 - Buffy coat (stored in 1 cryovial)
 - Whole blood (divided into 2 cryovials)
 - CSF (divided into 3 Cryovials)
 - Crypto isolate (in 2 Microbank vials)
 - qPCR whole blood (stored in 1 ceramic beaded tube)
 - qPCR plasma (stored in 1 eppendorf tube)
 - qPCR CSF supernatant (stored in 1 cryovial)
 - qPCR CSF pellet (stored in 1 cryovial)
 - qPCR sputum
 - qPCR urine
- Day 3:**
- Plasma (divided into 2 Cryovials)
 - Buffy coat (stored in 1 cryovial)
 - qPCR whole blood (stored in 1 ceramic beaded tube)
 - qPCR plasma (stored in 1 eppendorf tube)
 - qPCR urine
- Day 7:**
- Plasma (divided into 2 Cryovials)
 - Buffy coat (stored in 1 cryovial)
 - CSF (divided into 3 Cryovials)
 - Crypto isolate (in 2 Microbank vials – if growth from day 7 LP)
 - qPCR whole blood (stored in 1 ceramic beaded tube)
 - qPCR plasma (stored in 1 eppendorf tube)
 - qPCR CSF supernatant (stored in 1 cryovial)
 - qPCR CSF pellet (stored in 1 cryovial)
 - qPCR urine
- Day 14:**
- Plasma (divided into 2 Cryovials)
 - Buffy coat (stored in 1 cryovial)
 - CSF (divided in 2 Cryovials)
 - Crypto isolate (in 1 Microbank vial (if growth from Day 14 LP)
 - qPCR whole blood (stored in 1 ceramic beaded tube)
 - qPCR plasma (stored in 1 eppendorf tube)
 - qPCR CSF supernatant (stored in 1 cryovial)
 - qPCR CSF pellets (stored in 2 cryovial)

On D1, D3, D7 and D14 for qPCR sub-study please refer to the corresponding WPDs.

Lab Processing Chart – Cape Town

AMBITION-cm

High Dose AMBISOME on a Fluconazole Backbone for Cryptococcal Meningitis Induction Therapy in sub-Saharan Africa: A
Phase 3 Randomised Controlled Non-inferiority Trial

(As per Protocol version 2.1: 07/11/2017)

Protocol Version: 2.1 07/11/2017

Original LPC: Version 1.0, 14/06/2018

Revised/Updated: Version 1.2 13/09/2019: addition of qPCR sub-study in Cape Town

Laboratory Technologist: Kwana Lechiile

Email: klechiile@bhp.org.bw

Phone: +267 3902671

Lead Clinician: Dr David Lawrence

Clinical Advisor: Dr Timothee Boyer-Chammard

Laboratory Working Practice Document 3: Sample Processing and Storage

Table 1: Schedule of
(as per Study

| TABLE 6.1. Event Schedule | Screening | Week 1 | | | | | | | Week 2 | | | | | | | Wk 4 | Wk 6 | Wk 8 | Wk 10 | Wk 16 |
|-------------------------------------|-----------|--------|----|----|----|----|----|----|--------|----|-----|-----|-----|-----|-----|------|------|------|-------|-------|
| Study Day | ≤D0 | D1 | D2 | D3 | D4 | D5 | D6 | D7 | D8 | D9 | D10 | D11 | D12 | D13 | D14 | | | | | |
| Consent forms | | | | | | | | | | | | | | | | | | | | |
| PIS and signed consent | X | X | | | | | | | | | | | | | | | | | | |
| Follow up | | | | | | | | | | | | | | | | | | | | |
| Screening and randomisation | X | X | | | | | | | | | | | | | | | | | | |
| Clinical review | | X | X | X | X | X | X | X | X | X | X | X | X | X | X | | | | | |
| Outpatient follow-up | | | | | | | | | | | | | | | | X | X | X | X | |
| Week 16 telephone | | | | | | | | | | | | | | | | | | | | X |
| Clinical labs | | | | | | | | | | | | | | | | | | | | |
| HIV testing (if status unknown)* | | X | | | | | | | | | | | | | | | | | | |
| Pregnancy Test (Urine/Serum) (2) | | X | | | | | | | | | | | | | | | | | | |
| Full Blood Count | | X | | | | | | X | | | | | | | X | X | | | | |
| CD4 count & Viral load (if needed)* | | X | | | | | | | | | | | | | | | | | | |
| ALT | | X | | | | | | X | | | | | | | X | X | | | | |
| Urea, creatinine and electrolytes | | X | | X | | X | | X | | X | | X | | X | X | X | | | | |
| qPCR (4) | | X | | | | | | X | | | | | | | X | | | | | |
| Clinical evaluation | | | | | | | | | | | | | | | | | | | | |
| Chest X-ray (1) | | X | | | | | | | | | | | | | | | | | | |
| CSF | | | | | | | | | | | | | | | | | | | | |
| Opening pressure | | X | | | | | | X | | | | | | | X | | | | | |
| Cell count and differential* | | X | | | | | | | | | | | | | | | | | | |
| Protein, glucose* | | X | | | | | | | | | | | | | | | | | | |
| Routine culture * | | X | | | | | | | | | | | | | | | | | | |
| India ink examination* (3) | | X | | | | | | | | | | | | | | | | | | |
| Cryptococcal antigen* (3) | | X | | | | | | | | | | | | | | | | | | |
| Quantitative fungal culture | | X | | | | | | X | | | | | | | X | | | | | |
| Immune parameters and qPCR(4) | | X | | | | | | X | | | | | | | X | | | | | |


*Part of routine care. 1. If clinically indicated. 2. For women of childbearing age. 3. India ink or cryptococcal antigen required for inclusion. 4. As part of sub-studies at limited sites.

Events
Protocol)

Laboratory Working Practice Document 3: Sample Processing and Storage

Table 2: Specimen Log

| Evaluation | Specimen | Tests | Special Notes |
|----------------------|---------------------|--|---|
| Haematology | EDTA | FBC | All patients, when applicable |
| CD4 | EDTA | CD4 count | Should be performed if it has not been done within 3 months or if there is reason to suspect it is required. |
| Viral Load | EDTA | Viral load | Should be performed if it has not been done within 3 months or if there is reason to suspect it is required. |
| Resistance | EDTA | Resistance testing | If patient is on ART, has a detectable viral load, and reports good adherence or their ART history gives reason to suspect resistance, then resistance testing will be performed if possible at the site. |
| Chemistry | Yellow/Red Top Tube | ALT, Creatinine, Electrolytes | All patients, when applicable |
| Whole Blood | EDTA | For further analysis and sub-studies | Storage at -80°C until shipment |
| Plasma | EDTA | For further analysis and sub-studies | Storage at -80°C until shipment |
| Serum | Yellow/Red Top Tube | For further analysis and sub-studies | Storage at -80°C until shipment |
| Whole blood for qPCR | EDTA | For qPCR sub-study in blood | Storage at -80°C until shipment |
| Chemistry – CSF | CSF tube | Protein, Glucose | All patients, when applicable |
| Microbiology – CSF | CSF tube | India Ink, CrAg LFA, Gram stain, Cell count and differentiation, Routine culture (for D1 CSF and if clinically indicated thereafter) Quantitative cryptococcal culture (for all CSF) | Store isolate and any unused sample at -80°C until shipment |
| CSF for qPCR | CSF tube | For qPCR sub-study in CSF | Storage at -80°C until shipment |
| Sputum for qPCR | Sputum tube | For qPCR sub-study in sputum | Storage at -80°C until shipment |
| Urine for qPCR | Urine tube | For qPCR sub-study in urine | Storage at -80°C until shipment |

 **AMBITION-cm**
AMBIsome Therapy Induction OptimizatiON

Laboratory Working Practice Document 3: Sample Processing and Storage

Classification of samples:

Routine Tube: Part of routine clinical care, tests performed by usual laboratory and not required to be in an ISO accredited lab. Not curated by Ambition.

Monitoring Tube: As part of protocol driven monitoring, stored if samples remain and curated by Ambition.

Study Tube: For specified sub-studies or storage for future research. Curated by Ambition.

Laboratory Working Practice Document 3: Sample Processing and Storage

Table 3: Study lay-out

| Visit | Specimen | Tests Conducted | Storage | Special Notes | |
|---|---------------|--|--|---|--|
| <p>A number of samples on Day 1 are routine tests which may be handled by other laboratories *Store any unused samples* *Aim for a minimum of 5ml CSF from each LP*</p> | | | | | |
| Day 1 | CSF (5ml) x 3 | Routine Tube 1: perform protein and glucose, cell count and differential, CrAg LFA, India Ink, Gram stain, routine culture | No storage for study required. | May need to separate into two: biochemistry and microbiology | |
| | | Study Tube 1: perform quantitative cryptococcal culture in duplicates | Vortex and divide remaining CSF into 3x2ml cryovials (min. 200µl each) and store at -80°C | Store isolate once grown in 2 Microbank vials and keep at -80°C | |
| | | Study Tube 2: Sample for storage | General CSF storage: Unspun sample qPCR sub-study: centrifuge CSF at 2000g for 10 min and store supernatant and pellet in different cryovials at -80°C | Store at -80°C | |
| | EDTA (4ml) | Optional Routine Tube 1 and 2: CD4 and/or VL/Resistance Test (one tube for each test) | | | |
| | | Monitoring Tube 1: FBC | After FBC, spin at 800 ^x g for 15min and store remaining plasma as below | If being performed in research laboratory. | |
| | | Study Tube 1: Whole blood | Invert 8-10 times and pipette into 2x2ml cryovials (min. 900µl each) and store at -80°C | | |
| | | Study Tube 2: Spin tube at 800 ^x g for 15min and remove plasma | Divide plasma into 2x2ml cryovials. Remove buffy coat and put into 1 x 2ml cryovial | Store at -80°C | |
| | | Study Tube 3: qPCR sub-study | Remove 1.2ml and store in 2ml cryovial containing ceramic beads; centrifuge remaining 3.8ml at 2000g for 10min and store plasma in 1.5ml Eppendorf tube | Store at -80°C | |

Laboratory Working Practice Document 3: Sample Processing and Storage


| | | | | |
|-------|-------------------------------|---|---|--|
| | Yellow/Red Top Tube (4ml) x 2 | Monitoring Tube 1: ALT, urea, creatinine, and electrolytes | Spin and store remaining serum as below. | If being performed in research laboratory. |
| | | Study Tube 1: Spin tube at 800 ^x g for 15min and remove serum | Divide serum into 2x2ml cryovials and store at -80°C | |
| | Sputum | Sputum for qPCR sub-study | Process with dithiothreitol according to manufacturer's instructions (e.g. volume/volume and incubation 10min at 37°C). Centrifugation at 10,000 g for 3-5 min Remove the supernatant and discard Transfer the pellet into an 1-2ml Eppendorf tube | Store at -80°C |
| | Urine (20-50 ml) | Urine for qPCR sub-study | Centrifuge urine sample at 4000 rpm Remove the supernatant and discard Transfer the pellet into an 1-2ml Eppendorf tube | Store at -80°C |
| Day 3 | Yellow/Red Top Tube (4ml) | Monitoring Tube 1: Urea, creatinine and electrolytes | Spin and store remaining serum. | If being performed in research laboratory. |
| | EDTA (4ml) x 2 | Study Tube 1: Spin tube at 800 ^x g for 15min and remove plasma | Divide plasma into 2x2ml cryovials. Remove buffy coat and put into 1 x 2ml cryovial | Store at -80°C |
| | | Study Tube 2: qPCR sub-study | Remove 1.2ml and store in 2ml cryovial containing ceramic beads; centrifuge remaining 3.8ml at 2000g for 10min and store plasma in 1.5ml Eppendorf tube | Store at -80°C |
| | Urine (20-50 ml) | Urine for qPCR sub-study | Centrifuge urine sample at 4000 rpm Remove the supernatant and discard Transfer the pellet into an 1-2ml Eppendorf tube | Store at -80°C |
| Day 5 | Yellow/Red Top Tube (4ml) | Monitoring Tube 1: Urea, creatinine and electrolytes | Spin and store remaining serum. | If being performed in research laboratory. |

Laboratory Working Practice Document 3: Sample Processing and Storage

| | | | | |
|------------------|---------------------------|---|---|---|
| Day 7 | CSF (5ml) x 2 | Study Tube 1: Perform quantitative cryptococcal culture in duplicates | Vortex and divide remaining CSF into 3x2ml cryovials (min. 200µl each) and store at -80°C | Store isolate once grown in 2 Microbank vials and keep at -80°C |
| | | Study Tube 2: Sample for storage | General CSF storage: Unspun sample qPCR sub-study: centrifuge CSF at 2000g for 10min and store supernatant and pellet in different cryovials at -80°C | Store at -80°C |
| | EDTA (4ml) x 3 | Monitoring Tube 1: FBC | After FBC, spin at 800×g for 15min and store remaining plasma as below | If being performed in research laboratory. |
| | | Study Tube 1: Spin tube at 3500rpm/ 800×g for 15min and remove plasma | Divide plasma into 2x2ml cryovials. Remove buffy coat and put into 1 x 2ml cryovial | Store at -80°C |
| | | Study Tube 2: qPCR sub-study | Remove 1.2ml and store in 2ml cryovial containing ceramic beads; centrifuge remaining 3.8ml at 2000g for 10min and store plasma in 1.5ml Eppendorf tube | Store at -80°C |
| | Yellow/Red Top Tube (4ml) | Monitoring Tube 1: ALT, urea, creatinine, and electrolytes | Spin and store remaining serum. | If being performed in research laboratory. |
| Urine (20-50 ml) | Urine for qPCR sub-study | Centrifuge urine sample at 4000 rpm Remove the supernatant and discard Transfer the pellet into an 1-2ml Eppendorf tube | Store at -80°C | |
| Day 10 | Yellow/Red Top Tube (4ml) | Urea, creatinine and electrolytes | Spin and store remaining serum. | If being performed in research laboratory. |
| Day 12 | Yellow/Red Top Tube (4ml) | Monitoring Tube 1: Urea, creatinine and electrolytes | Spin and store remaining serum. | If being performed in research laboratory. |
| Day 14 | CSF (5ml) x 2 | Study Tube 1: Perform quantitative cryptococcal culture in duplicates | Vortex and divide remaining CSF into 3x2ml cryovials (min. 200µl each) and store at -80°C | Store isolate once grown in 2 Microbank vials and keep at -80°C |
| | | Study Tube 2: Sample for storage | General CSF storage: | Store at -80°C |

Laboratory Working Practice Document 3: Sample Processing and Storage

| | | | | |
|--------|---------------------------|--|---|--|
| | | | <p>Unspun sample</p> <p>India ink storage: centrifuge 500µl of CSF at 2000g for 10 min and store the pellet at -80°C.</p> <p>qPCR sub-study: centrifuge all remaining CSF at 2000g for 10min and store supernatant and pellet in different cryovials at -80°C</p> | |
| | EDTA (4ml) x 3 | Monitoring Tube 1: FBC | After FBC, spin at 800×g for 15min and store remaining plasma as below | If being performed in research laboratory. |
| | | Study Tube 1: Spin tube at 800×g for 15min and remove plasma | Divide plasma into 2x2ml cryovials. Remove buffy coat and put into 1 x 2ml cryovial | Store at -80°C |
| | | Study Tube 2: qPCR sub-study | Remove 1.2ml and store in 2ml cryovial containing ceramic beads; centrifuge remaining 3.8ml at 2000g for 10min and store plasma in 1.5ml Eppendorf tube | Store at -80°C |
| | Yellow/Red Top Tube (4ml) | Monitoring Tube 1: ALT, urea, creatinine, and electrolytes | Spin and store remaining serum. | If being performed in research laboratory. |
| Day 28 | EDTA (4ml) | Monitoring Tube 1: FBC | After FBC, spin at 800×g for 15min and store remaining plasma at -80°C | If being performed in research laboratory. |
| | Yellow/Red Top Tube (4ml) | Monitoring Tube 1: ALT, urea, creatinine, and electrolytes | Spin and store remaining serum at -80°C. | If being performed in research laboratory. |



Laboratory Working Practice Document 3: Sample Processing and Storage

Table 4: Shipments

| Specimen | Destination | Shipment Date | Special Notes |
|--------------------|-------------|---------------|---------------|
| CSF | | | |
| CSF – PK/PD Study | | | |
| Crypto Isolates | | | |
| CSF and Blood qPCR | | | |
| PK Plasma | | | |
| Serum | | | |
| Whole Blood | | | |
| Plasma | | | |
| Buffy Coat | | | |
| Sputum qPCR | | | |
| Urine qPCR | | | |

Laboratory Working Practice Document 3: Sample Processing and Storage

Training

Each staff member receives or has direct access to applicable Working Practice Documents (WPDs).

Each staff member reviews the applicable WPDs once a year.

All WPD training is documented and tracked in the training log located in the Investigator Site File (ISF)

New staff is trained on applicable WPDs within 30 days of employment and all WPDs within 90 days of employment.

Staff members whose duties fall within this WPD scope are retrained within 14 days of the approval of each WPD revision.

Laboratory Working Practice Document 3: Sample Processing and Storage



Staff signatures: (signing below indicates that you have read this SOP and understand the material contained in it)

| Date | Name (Please print) | Signature |
|------|---------------------|-----------|
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |
| | | |