





Laboratory Working Practice Document 3: Sample Processing and Storage

Ambition Trial Coordinating Centre
Private Bag 320, Princess Marina Hospital
Gaborone, Botswana

Laboratory Work Practice Document: 3 (Gaborone) Sample Processing and Storage

Title of study	High Dose AMBISOME® on a Fluconazole Backbone for Cryptococcal Meningitis Induction Therapy in sub-Saharan Africa: A Phase III Randomized Controlled Non-inferiority Trial		
Acronym	Ambition-cm – AMBIsome Therapy Induction Optimization		
ISRCTN No.:	ISRCTN72509687		
WPD Current version	Version 1.3 22/06/2020		
Author(s)	Kwana Lechiile Lab Scientist		22/06/2020
	Timothée Boyer Chamard Clinical Advisor		22/06/2020
Reviewer(s)	David Lawrence Lead Clinician		22/06/2020
Approved by	Joseph Jarvis CI		22/06/2020

Revision History:		
Version Number	Effective Date	Reason for Change
1.0	14/06/2018	First version
1.1	22/08/2018	Addition of India ink storage at D14
1.2	10/09/2019	Addition of PAXgene Specimen Collection and Storage
1.3	22/06/2020	Addition of PBMC isolation

Laboratory Working Practice Document 3: Sample Processing and Storage

Purpose

This document describes the day to day laboratory duties

References

Ambition Phase III Trial Protocol

Materials

Lab WPD 1 – Laboratory equipment
Lab WPD 2 – Quantitative Cryptococcal cultures
Lab WPD 7 – Quantitative Cryptococcal PCR sub-study
Lab WPD 8 – Cryptococcal Semi-Quantitative Antigen
Lab WPD 9 - PAXgene Specimen Collection and Storage
Clinical WPD 15- Timing of evaluation and tests
Clinical WPD 20 – Quantitative PCR Sub-study

Appendices

Lab Processing Chart: Gaborone

Section A – Lab maintenance, quality control

Monitoring of freezer and fridge temperatures, reagents/agar expiration dates, working order of autoclaves, incubators and other laboratory equipment, should all be performed in line with on-site laboratory quality control and health and safety procedures

Weekly stock check of lab consumables for trial. Inform trial team on-site and trial manager when replacements required. See Laboratory WPD 1 – Laboratory equipment

Laboratory Working Practice Document 3: Sample Processing and Storage

Section B – Samples required

All results should be recorded with patient study number, date, and study day.

All results should be communicated to the study team.

– Study bloods

The following must be performed according to the trial protocol

- CD4 at baseline if not done within 3 months
- FBC and ALT on days 1, 7, 14, 28
- Creatinine, electrolytes on days 1, 3, 5, 7, 10, 12, 14, 28
- Biosynex CryptoPS CrAg on whole blood on D1
(please see Lab WPD 8 – Cryptococcal Semi-Quantitative Antigen)
- qPCR samples (4ml of blood – EDTA tubes): D1 (pre-dose), 24 hours post-dose, D3, D7 and D14,
(please see Appendix, Lab WPD7 & Clinical WPD 20 - Quantitative PCR Sub-study)
- PAXgene tube on D1 (2.5mL of blood)
(please see Appendix and Lab WPD 9 - PAXgene Specimen Collection and Storage)
- PBMC isolation on D1 (20mL of blood)

– CSF

- India Ink examination and Cryptococcal antigen (IMMY CrAg) on D1
- Biosynex CryptoPS CrAg on CSF on D1
(please see Lab WPD 8 – Cryptococcal Semi-Quantitative Antigen)
- CSF cell count and differential, protein, glucose and routine culture on D1 for all patients and additional LPs, as indicated
- Quantitative fungal culture will be conducted on D1, 7, and 14
- CSF for qPCR (1ml): D1 (pre-dose), D7 and D14
(please see Appendix, Lab WPD7 & Clinical WPD 20 - Quantitative PCR Sub-study)
- PAXgene tube on D1 (2.5mL of CSF)
(please see Appendix and Lab WPD 9 - PAXgene Specimen Collection and Storage)

Laboratory Working Practice Document 3: Sample Processing and Storage

Table 1: Timings of study bloods and CSF analyses

Event Schedule	Screening	Week 1							Week 2							Wk 4
Study Day	≤D0	D1	D2	D3	D4	D5	D6	D7	D8	D9	D10	D11	D12	D13	D14	
Clinical labs																
HIV testing (If status unknown)*	X															
Pregnancy Test (Urine/Serum) (1)	X															
Full Blood Count	X							X							X	X
CD4 count & Viral load (If needed)*	X															
ALT	X							X							X	X
Urea, creatinine and electrolytes	X		X		X		X			X		X		X	X	X
Blood for drug levels (2)	X								X							
CSF																
Opening pressure	X							X							X	
Cell count and differential*	X															
Protein, glucose*	X															
Routine culture *	X															
India Ink examination* (3)	X															
Cryptococcal antigen* (3)	X															
Quantitative fungal culture	X							X							X	
CSF Drug levels (2)	X							X							X	
Immune parameters (2)	X							X							X	

*Part of routine care. 1. For women of childbearing age. 2. Sub-studies at limited sites 3. India Ink or CrAg required for inclusion.

Laboratory Working Practice Document 3: Sample Processing and Storage

Section C - Procedures for sample processing and saving - Please see Appendix: Lab Processing Chart

Day 1 (or 0)

Required samples:

- 2 yellow (or red) top (1 for biochemistry, 1 for serum storage)
- ≥3 purple top (1 for FBC, 1 for whole blood storage, 1 for plasma and buffy coat storage and other samples may be for CD4, viral load and resistance testing)
- CSF sample in 3 x white top tube (5 ml each)
- qPCR samples
(Please see WPD 20 – Quantitative PCR sub-study)
- 2 PAXgene tubes (1 for blood, 1 for CSF)
(Please see WPD 9 - PAXgene Specimen Collection and Storage)
- 20mL lithium heparin blood for PBMC isolation

DAY 1 CSF SAMPLES

NB CSF must be processed within **4 HOURS** of receipt

1. Perform protein and glucose on CSF
2. Perform India Ink, CrAg (if required), cell count (on day 0/1 sample)
3. Plate out day 1 unspun CSF dilutions 0-4 (neat, one in 10, 100, 1,000 and 10,000) for quantitative culture (see WPD 2)

Store isolate once grown in 2 Microbank vials.
4. Vortex CSF and save: divide into 3 x 2ml Cryovials - min 200 µL each. Label and place into freezer box.
5. qPCR SAMPLES: Please see Lab WPD7 & Clinical WPD 20 – Quantitative PCR Sub-study
6. Biosynex CryptoPS CrAg: Please see Lab WPD 8 – Cryptococcal semi-quantitative antigen
7. PAXgene tube: Please see WPD 9 - PAXgene Specimen Collection and Storage

TB culture only if clinically indicated (will be performed in TB lab, and is not mandated by trial protocol)

DAY 1 BLOOD SAMPLES

1. Take blood for CD4, viral load and/or resistance testing (if indicated), whole blood, plasma and buffy coat (all purple tops) and serum (yellow or red top). Total required = 2 yellow or red top (4 ml) and ≥3 purple tops (4 ml each)
2. Perform FBC, ALT, Creatinine and electrolytes on Day 1
3. **WHOLE BLOOD:** take 2x cryovials and pipette 900 uL of WHOLE BLOOD (unspun sample) into each (should be mixed well in purple top tubes prior to transfer). Label and place into whole blood storage box.

Laboratory Working Practice Document 3: Sample Processing and Storage

4. PLASMA AND BUFFY COAT: Spin blood from purple top tube at 800xg for 15 minutes. Avoiding the buffy coat and red cells at the bottom, remove as much of the plasma as possible using a sterile Pasteur pipette and divide into 2 cryovials.

Buffy coat is the small region between the plasma and red blood cells. After removing the plasma, use a sterile Pasteur pipette to carefully aspirate the buffy coat in a circular movement trying not to disturb the red blood cell layer then store in a single cryovial

5. SERUM: Spin yellow (or red) top sample at 800xg for 15 minutes. In the yellow top, gel now separates cells from serum. . Use a sterile Pasteur pipette to remove as much of the serum possible and divide into 2 cryovials.
6. qPCR SAMPLES: Please see Lab WPD7 & Clinical WPD 20 – Quantitative PCR Sub-study
7. Biosynex CryptoPS CrAg: Please see Lab WPD 8 – Cryptococcal semi-quantitative antigen
8. PAXgene tube: Please see WPD 9 - PAXgene Specimen Collection and Storage
9. PBMC Isolation: 20mL blood in Lithium Heparin tube

24 hours post-dose

Required samples:

- qPCR blood: Please see WPD 20 – Quantitative PCR Sub-study

Day 3

Required samples:

- 1 purple top (for plasma and buffy coat)
 - 1 yellow/red top tube (4ml)
 - qPCR samples
- (Please see Lab WPD7 & Clinical WPD 20 – Quantitative PCR sub-study)

DAY 3 BLOOD SAMPLES

1. Perform urea, creatinine and electrolytes on Day 3
2. PLASMA AND BUFFY COAT
3. qPCR SAMPLES: Please see Lab WPD7 & Clinical WPD 20 – Quantitative PCR Sub-study

Laboratory Working Practice Document 3: Sample Processing and Storage

Day 5

Required samples:

- 1 yellow/red top tube (4ml)

DAY 5 BLOOD SAMPLES

1. Perform urea, creatinine and electrolytes on Day 5

Day 7

Required samples:

- 2 purple top (1 for FBC, 1 for plasma and buffy coat storage)
- 1 yellow/red top tube (4ml)
- CSF sample in 2 x white top tubes (5 ml each)
- qPCR samples (Please see Lab WPD7 & Clinical WPD 20 – Quantitative PCR sub-study)

DAY 7 CSF SAMPLES

1. Perform quantitative cultures on Day 7 LP.
2. Save remaining CSF, divide into 3 x 2ml Cryovials.
3. Save Cryptococcus isolate in 2 Microbank vials if growth from Day 7 LP
4. qPCR SAMPLES: Please see Lab WPD7 & Clinical WPD 20 – Quantitative PCR Sub-study

DAY 7 BLOOD SAMPLES

1. Perform FBC, ALT, urea, creatinine and electrolytes on Day 7
2. PLASMA AND BUFFY COAT
3. qPCR SAMPLES: Please see Lab WPD7 & Clinical WPD 20 – Quantitative PCR Sub-study

Day 10

Required samples:

- 1 yellow/red top tube (4ml)

DAY 10 BLOOD SAMPLES

1. Perform urea, creatinine and electrolytes on Day 10

Laboratory Working Practice Document 3: Sample Processing and Storage

Day 12

Required samples:

- 1 yellow/red top tube (4ml)

DAY 12 BLOOD SAMPLES

1. Perform urea, creatinine and electrolytes on Day 12

Day 14

Required samples:

- 2 purple top (1 for FBC and 1 for plasma and buffy coat)
 - 1 yellow/red top tube (4ml)
 - CSF sample in 2 x white top tube (5 mls each)
 - qPCR samples (including India ink storage)
- (Please see Lab WPD7 & Clinical WPD 20 – Quantitative PCR sub-study)

DAY 14 CSF SAMPLES

1. Perform quantitative cultures on Day 14 LP. Save remaining CSF, divide into 3 x 2ml Cryovials.
2. Save Cryptococcus isolate in 2 Microbank vial if growth from Day 14 LP
3. qPCR SAMPLES: Please see WPD 20 – Quantitative PCR Sub-study

DAY 14 BLOOD SAMPLES

1. Perform FBC, ALT, urea, creatinine and electrolytes on Day 14
2. PLASMA AND BUFFY COAT
3. qPCR SAMPLES: Please see WPD 20 – Quantitative PCR Sub-study

Day 28

Required samples:

- 1 purple top (for FBC)
- 1 yellow/red top tube (4ml)

DAY 28 BLOOD SAMPLES

1. Perform FBC, ALT, urea, creatinine and electrolytes on Day 28

Laboratory Working Practice Document 3: Sample Processing and Storage

PROCESSING AND LABELLING

SPINNING: For all samples use 800xg for 15 mins, unless stated otherwise on LPC.

LABELLING: Each sample must be labelled with study number (e.g. AMB 1-1-001), nature of sample (e.g. CSF/plasma etc), study day (e.g. d3), date and time (time specimen taken, not time specimen received in lab).

The study day is derived from **the first day the patient receives study drugs**. That is day 1. That is not necessarily the day of hospital admission.

STORAGE: Cryovials containing samples should be stored in appropriately labelled cryoboxes. Cryoboxes should be stored in -80°C freezer

CRYPTOCOCCAL ISOLATE STORAGE: Cryptococcal isolates are stored on Microbank beads in an appropriately labelled cryobox and stored in -80°C freezer.

PROCEDURE FOR STORING ON MICROBANK BEADS:

1. Using a permanent marker, label with the study number, date and day of study (e.g. AMB 1-1-001, 12.10.15, d7)
2. Under aseptic conditions open the screw cap cryovial.
3. Inoculate the cryopreservative fluid with young colonial growth (18-24 hours) picked from a pure culture to approximately a 3-4 McFarland standard.
4. Close vial tightly and invert 4-5 times to emulsify organism. **DO NOT VORTEX!**
5. At this point the microorganisms will be bound to the porous beads. The excess cryopreservative fluid should be well aspirated leaving the inoculated beads as free of liquid as possible. Close the vial finger tight.
6. Record on the storage log.
7. Store the inoculated cryovial in -80°C freezer.

Laboratory Working Practice Document 3: Sample Processing and Storage

Section D - Summary of samples to be saved

- Day 1:** Serum (stored in 2 Cryovials)
Plasma (divided into 2 Cryovials)
Buffy coat (stored in 1 cryovial)
Whole blood (divided into 2 cryovials)
CSF (divided into 3 Cryovials)
Crypto isolate (in 2 Microbank vials)
qPCR whole blood (stored in 1 ceramic beaded tube)
qPCR plasma (stored in 1 eppendorf tube)
qPCR CSF supernatant (stored in 1 cryovial)
qPCR CSF pellet (stored in 1 cryovial)
Blood PAXgene tube
CSF PAXgene tube
PBMCs
- 24HR:** qPCR whole blood (stored in 1 ceramic beaded tube)
qPCR plasma (stored in 1 eppendorf tube)
- Day 3:** Plasma (divided into 2 Cryovials)
Buffy coat (stored in 1 cryovial)
qPCR whole blood (stored in 1 ceramic beaded tube)
qPCR plasma (stored in 1 eppendorf tube)
- Day 7:** Plasma (divided into 2 Cryovials)
Buffy coat (stored in 1 cryovial)
CSF (divided into 3 Cryovials)
Crypto isolate (in 2 Microbank vials – if growth from day 7 LP)
qPCR whole blood (stored in 1 ceramic beaded tube)
qPCR plasma (stored in 1 eppendorf tube)
qPCR CSF supernatant (stored in 1 cryovial)
qPCR CSF pellet (stored in 1 cryovial)
- Day 14:** Plasma (divided into 2 Cryovials)
Buffy coat (stored in 1 cryovial)
CSF (divided in 2 Cryovials)
Crypto isolate (in 1 Microbank vial (if growth from Day 14 LP)
qPCR whole blood (stored in 1 ceramic beaded tube)
qPCR plasma (stored in 1 eppendorf tube)
qPCR CSF supernatant (stored in 1 cryovial)
qPCR CSF pellets (stored in 2 cryovial)
- On D1, D3, D7 and D14 for qPCR or semi-quantitative CrAg sub-studies please refer to the corresponding WPDs.**

Lab Processing Chart - Botswana

AMBITION-cm

High Dose AMBISOME on a Fluconazole Backbone for Cryptococcal Meningitis Induction Therapy in sub-Saharan Africa: A
Phase 3 Randomised Controlled Non-inferiority Trial

(As per Protocol version 2.1: 07/11/2017)

Protocol Version: 2.1 07/11/2017

Original LPC: Version 1.0, 14/06/2018

Revised/Updated: Version 1.1 03/07/2018: addition of India ink storage at D14

Version 1.2 10/09/2019 : addition of PAXgene specimen collection and storage at D1

Laboratory Technologist: Kwana Lechiile

Email: klechiile@bhp.org.bw

Phone: +267 3902671

Lead Clinician: Dr David Lawrence

Clinical Advisor: Dr Timothee Boyer-Chammard

Laboratory Working Practice Document 3: Sample Processing and Storage

Table 1: Schedule of Events
(as per Study Protocol)


TABLE 6.1. Event Schedule	Screening	Week 1							Week 2							Wk 4	Wk 6	Wk 8	Wk 10	Wk 16
Study Day	≤D0	D1	D2	D3	D4	D5	D6	D7	D8	D9	D10	D11	D12	D13	D14					
Consent forms																				
PIS and signed consent	X	X																		
Follow up																				
Screening and randomisation	X	X																		
Clinical review		X	X	X	X	X	X	X	X	X	X	X	X	X	X					
Outpatient follow-up																X	X	X	X	
Week 16 telephone																				X
Clinical labs																				
HIV testing (if status unknown)*		X																		
Pregnancy Test (Urine/Serum) (2)		X																		
Full Blood Count		X						X						X		X				
CD4 count & Viral load (if needed)*		X																		
ALT		X						X						X		X				
Urea, creatinine and electrolytes		X		X		X		X		X		X		X		X				
qPCR (4)		X						X						X						
Clinical evaluation																				
Chest X-ray (1)		X																		
CSF																				
Opening pressure		X						X						X						
Cell count and differential*		X																		
Protein, glucose*		X																		
Routine culture *		X																		
India ink examination* (3)		X																		
Cryptococcal antigen* (3)		X																		
Quantitative fungal culture		X						X						X						
Immune parameters and qPCR(4)		X						X						X						

*Part of routine care. 1. If clinically indicated. 2. For women of childbearing age. 3. India ink or cryptococcal antigen required for inclusion. 4. As part of sub-studies at limited sites.

Laboratory Working Practice Document 3: Sample Processing and Storage

Table 2: Specimen Log

Evaluation	Specimen	Tests	Special Notes
Haematology	EDTA	FBC	All patients, when applicable
CD4	EDTA	CD4 count	Should be performed if it has not been done within 3 months or if there is reason to suspect it is required.
Viral Load	EDTA	Viral load	Should be performed if it has not been done within 3 months or if there is reason to suspect it is required.
Resistance	EDTA	Resistance testing	If patient is on ART, has a detectable viral load, and reports good adherence or their ART history gives reason to suspect resistance, then resistance testing will be performed if possible at the site.
Chemistry	Yellow/Red Top Tube	ALT, Creatinine, Electrolytes	All patients, when applicable
Whole Blood	EDTA	For further analysis and sub-studies	Storage at -80°C until shipment
Plasma	EDTA	For further analysis and sub-studies	Storage at -80°C until shipment
Serum	Yellow/Red Top Tube	For further analysis and sub-studies	Storage at -80°C until shipment
Whole blood for qPCR	EDTA	For qPCR sub-study in blood	Storage at -80°C until shipment
Chemistry – CSF	CSF tube	Protein, Glucose	All patients, when applicable
Microbiology – CSF	CSF tube	India Ink, CrAg LFA, Gram stain, Cell count and differentiation, Routine culture (for D1 CSF and if clinically indicated thereafter) Quantitative cryptococcal culture (for all CSF)	Store isolate and any unused sample at -80°C until shipment
CSF for qPCR	CSF tube	For qPCR sub-study in CSF	Storage at -80°C until shipment
Blood PAXgene	PAXgene tube	For PAXgene tube sub-study	Storage immediately at 18-25°C, then at -20°C or at -80°C
CSF PAXgene	PAXgene tube	For PAXgene tube sub-study	Storage immediately at 18-25°C, then at -20°C or at -80°C
PBMCs	Lithium heparin	For PBMC isolation	As per local SOP.

 **AMBITION-cm**
AMBIsome Therapy Induction OptimizatiON

Laboratory Working Practice Document 3: Sample Processing and Storage

Classification of samples:

Routine Tube: Part of routine clinical care, tests performed by usual laboratory and not required to be in an ISO accredited lab. Not curated by Ambition.

Monitoring Tube: As part of protocol driven monitoring, stored if samples remain and curated by Ambition.

Study Tube: For specified sub-studies or storage for future research. Curated by Ambition.

Laboratory Working Practice Document 3: Sample Processing and Storage

Table 3: Study lay-out

Visit	Specimen	Tests Conducted	Storage	Special Notes
<p>A number of samples on Day 1 are routine tests which may be handled by other laboratories *Store any unused samples* *Aim for a <u>minimum</u> of 5ml CSF from each LP*</p>				
Day 1	CSF (5ml) x 3	Routine Tube 1: perform protein and glucose, cell count and differential, CrAg LFA, India Ink, Gram stain, routine culture	No storage for study required.	May need to separate into two: biochemistry and microbiology
		Study Tube 1: perform quantitative cryptococcal culture in duplicates	Vortex and divide remaining CSF into 3x2ml cryovials (min. 200µl each) and store at -80°C	Store isolate once grown in 2 Microbank vials and keep at -80°C
		Study Tube 2: Sample for storage	General CSF storage: Unspun sample qPCR sub-study: centrifuge CSF at 2000g for 10 min and store supernatant and pellet in different cryovials at -80°C	Store at -80°C
	CSF PAXgene tube	Immediately after adding CSF, gently invert the PAXgene tube 8-10 times. Samples do not require any processing.	Immediately after collecting specimens, store the PAXgene tube at 18-25°C for a minimum of 2 hours and a maximum of 72 hours. The PAXgene tubes can then be stored at -20°C or below. If you would prefer to store at -80°C then it will be necessary to first store at -20°C for at least 24 hours before transferring to the deep freeze.	When storing in the freezer the preference is to store on a wire rack: do not freeze tubes upright in a Styrofoam tray as this may cause the tubes to crack. If you cannot store on a wire rack then keep the tubes in a freezer bag and store them in a safe, secure position in the freezer.
	EDTA (4ml)	Optional Routine Tube 1 and 2: CD4 and/or VL/Resistance Test (one tube for each test)		
Monitoring Tube 1: FBC			After FBC, spin at 800×g for 15min and store remaining plasma as below	If being performed in research laboratory.

Laboratory Working Practice Document 3: Sample Processing and Storage


		Study Tube 1: Whole blood	Invert 8-10 times and pipette into 2x2ml cryovials (min. 900µl each) and store at -80°C	
		Study Tube 2: Spin tube at 800*g for 15min and remove plasma	Divide plasma into 2x2ml cryovials. Remove buffy coat and put into 1 x 2ml cryovial	Store at -80°C
		Study Tube 3: qPCR sub-study	Remove 1.2ml and store in 2ml cryovial containing ceramic beads; centrifuge remaining 3.8ml at 2000g for 10min and store plasma in 1.5ml Eppendorf tube	
	Yellow/Red Top Tube (4ml) x 2	Monitoring Tube 1: ALT, urea, creatinine, and electrolytes	Spin and store remaining serum as below.	If being performed in research laboratory.
		Study Tube 1: Spin tube at 800*g for 15min and remove serum	Divide serum into 2x2ml cryovials and store at -80°C	
	Blood PAXgene tube	Immediately after adding the blood, gently invert the PAXgene tube 8-10 times. Samples do not require any processing.	Immediately after collecting specimens, store the PAXgene tube at 18-25°C for a minimum of 2 hours and a maximum of 72 hours. The PAXgene tubes can then be stored at -20°C or below. If you would prefer to store at -80°C then it will be necessary to first store at -20°C for at least 24 hours before transferring to the deep freeze.	When storing in the freezer the preference is to store on a wire rack: do not freeze tubes upright in a Styrofoam tray as this may cause the tubes to crack. If you cannot store on a wire rack then keep the tubes in a freezer bag and store them in a safe, secure position in the freezer.
	Lithium heparin (20ml)	PBMC Isolation	As per local SOP	
24 hours post-dose	EDTA (4ml) x 1	Study Tube 1: qPCR sub-study	Remove 1.2ml and store in 2ml cryovial containing ceramic beads; centrifuge remaining 3.8ml at 2000g for 10min and store plasma in 1.5ml Eppendorf tube	qPCR tube is collected within 1 h of t23
Day 3	Yellow/Red Top Tube (4ml)	Monitoring Tube 1: Urea, creatinine and electrolytes	Spin and store remaining serum.	If being performed in research laboratory.

Laboratory Working Practice Document 3: Sample Processing and Storage

	EDTA (4ml) x 2	Study Tube 1: Spin tube at 800*g for 15min and remove plasma	Divide plasma into 2x2ml cryovials. Remove buffy coat and put into 1 x 2ml cryovial	Store at -80°C
		Study Tube 2: qPCR sub-study	Remove 1.2ml and store in 2ml cryovial containing ceramic beads; centrifuge remaining 3.8ml at 2000g for 10min and store plasma in 1.5ml Eppendorf tube	
Day 5	Yellow/Red Top Tube (4ml)	Monitoring Tube 1: Urea, creatinine and electrolytes	Spin and store remaining serum.	If being performed in research laboratory.
Day 7	CSF (5ml) x 2	Study Tube 1: Perform quantitative cryptococcal culture in duplicates	Vortex and divide remaining CSF into 3x2ml cryovials (min. 200µl each) and store at -80°C	Store isolate once grown in 2 Microbank vials and keep at -80°C
		Study Tube 2: Sample for storage	General CSF storage: Unspun sample qPCR sub-study: centrifuge CSF at 2000g for 10min and store supernatant and pellet in different cryovials at -80°C	Store at -80°C
	EDTA (4ml) x 3	Monitoring Tube 1: FBC	After FBC, spin at 800*g for 15min and store remaining plasma as below	If being performed in research laboratory.
		Study Tube 1: Spin tube at 3500rpm/ 800*g for 15min and remove plasma	Divide plasma into 2x2ml cryovials. Remove buffy coat and put into 1 x 2ml cryovial	Store at -80°C
		Study Tube 2: qPCR sub-study	Remove 1.2ml and store in 2ml cryovial containing ceramic beads; centrifuge remaining 3.8ml at 2000g for 10min and store plasma in 1.5ml Eppendorf tube	
	Yellow/Red Top Tube (4ml)	Monitoring Tube 1: ALT, urea, creatinine, and electrolytes	Spin and store remaining serum.	If being performed in research laboratory.
Day 10	Yellow/Red Top Tube (4ml)	Urea, creatinine and electrolytes	Spin and store remaining serum.	If being performed in research laboratory.

Laboratory Working Practice Document 3: Sample Processing and Storage

Day 12	Yellow/Red Top Tube (4ml)	Monitoring Tube 1: Urea, creatinine and electrolytes	Spin and store remaining serum.	If being performed in research laboratory.
Day 14	CSF (5ml) x 2	Study Tube 1: Perform quantitative cryptococcal culture in duplicates	Vortex and divide remaining CSF into 3x2ml cryovials (min. 200µl each) and store at -80°C	Store isolate once grown in 2 Microbank vials and keep at -80°C
		Study Tube 2: Sample for storage	General CSF storage: Unspun sample India ink storage: centrifuge 500µl of CSF at 2000g for 10 min and store the pellet at -80°C. qPCR sub-study: centrifuge all remaining CSF at 2000g for 10min and store supernatant and pellet in different cryovials at -80°C	Store at -80°C
	EDTA (4ml) x 3	Monitoring Tube 1: FBC	After FBC, spin at 800×g for 15min and store remaining plasma as below	If being performed in research laboratory.
		Study Tube 1: Spin tube at 800×g for 15min and remove plasma	Divide plasma into 2x2ml cryovials. Remove buffy coat and put into 1 x 2ml cryovial	Store at -80°C
		Study Tube 2: qPCR sub-study	Remove 1.2ml and store in 2ml cryovial containing ceramic beads; centrifuge remaining 3.8ml at 2000g for 10min and store plasma in 1.5ml Eppendorf tube	
	Yellow/Red Top Tube (4ml)	Monitoring Tube 1: ALT, urea, creatinine, and electrolytes	Spin and store remaining serum.	If being performed in research laboratory.
Day 28	EDTA (4ml)	Monitoring Tube 1: FBC	After FBC, spin at 800×g for 15min and store remaining plasma at -80°C	If being performed in research laboratory.
	Yellow/Red Top Tube (4ml)	Monitoring Tube 1: ALT, urea, creatinine, and electrolytes	Spin and store remaining serum at -80°C.	If being performed in research laboratory.



Laboratory Working Practice Document 3:
Sample Processing and Storage

Table 4: Shipments

Specimen	Destination	Shipment Date	Special Notes
CSF			
CSF – PK/PD Study			
Crypto Isolates			
CSF and Blood qPCR			
PK Plasma			
Serum			
Whole Blood			
Plasma			
Buffy Coat			
CSF and Blood PAXgene tubes			
PBMCs			

Laboratory Working Practice Document 3: Sample Processing and Storage

Training

Each staff member receives or has direct access to applicable Working Practice Documents (WPDs).

Each staff member reviews the applicable WPDs once a year.

All WPD training is documented and tracked in the training log located in the Investigator Site File (ISF)

New staff is trained on applicable WPDs within 30 days of employment and all WPDs within 90 days of employment.

Staff members whose duties fall within this WPD scope are retrained within 14 days of the approval of each WPD revision.

Laboratory Working Practice Document 3: Sample Processing and Storage

Staff signatures: (signing below indicates that you have read this SOP and understand the material contained in it)

Date	Name (Please print)	Signature